

## Lab: Preparation of Bromoethane and Investigation of its Properties

### Reaction



In this preparation you will meet the techniques of reflux, simple distillation and use of a separating funnel, which are commonly used in the purification of an organic liquid to give the best possible yield.

### Safety

Ethanol is flammable. Concentrated sulfuric acid is corrosive. Therefore, you must: keep stoppers on bottles as much as possible, keep flammable liquids away from flames, and wear safety glasses.

### Procedure

1. In the round bottomed flask, place 7 cm<sup>3</sup> of ethanol, 3 cm<sup>3</sup> of conc. H<sub>2</sub>SO<sub>4</sub> and 12 g of KBr. Add a boiling stone and reflux gently for 10 minutes.
2. Set up the apparatus for distillation and distill off the bromoethane, collecting it under water in a small beaker (bp 38°C).

### See Diagrams

**Refluxing:** is the process of boiling volatile reactants in a flask connected to a condenser which is fitted vertically so that the condensed liquid runs back into the flask. In this way, the liquid mixture can be maintained at a fairly high temperature for a long period of time without loss of evaporation.

### Tests

1. Hydrolysis
  - Boil 2-3 drops with a little dilute nitric acid then add aqueous AgNO<sub>3</sub> to show the formation of a halogen.
2. Reaction with alcoholic potash
  - In a boiling tube fitted with a cork and bent delivery tube place 10 drops of liquid, 5 cm<sup>3</sup> of alcohol and 5-6 pellets of KOH. Warm gently and collect 3 test tubes of gas over water. Test the gas...
    - a) by burning
    - b) with bromine water (in TTE)
    - c) with acidified KMnO<sub>4</sub>

What forms each time?

Write equations to represent the reactions occurring.