Oxidation States of Manganese

Introduction

Atoms in compounds may be found in various oxidation states. Atoms of some elements, like manganese, have many possible oxidation states. Oxidation-reduction reactions lead to changes in oxidation states.

Purpose

The oxidation state of manganese in potassium permanganate is changed through a series of redox reactions. The redox reactions are controlled by adjusting the acidity of the reaction medium.

Safety

- 1. Wear protective goggles throughout the laboratory activity.
- 2. Dilute sulfuric acid and dilute sodium hydroxide solutions are corrosive. Solid potassium permanganate is a strong oxidant. Potassium permanganate and sodium bisulfite, (aka: sodium hydrogen sulphite, NaHSO₃), are toxic.
- 3. Use caution when handling the dilute sulfuric acid solution.
- 4. Do not permit the potassium permanganate to come in contact with any reducing agents. Save any spilled material for disposal according to the procedure given.
- 5. Dispose of all materials as your teacher directs.

Procedure

Use pulled Beral pipets.

- 1. Put 10 drops KMnO₄ solution, (the oxidising agent, why?), in four wells of a 24-well plate, A-1 through A-4. (The well A-1 is your control)
- 2. Add four drops NaOH solution to well A-2, thus basic medium.
- 3. Add three drops 3M H₂SO₄ solution top well A-4, identifying an acidic medium.
- 4. Place the plate on a sheet of white paper or other white background, for visibility.
- 5. Add three drops NaHSO₃ solution, (the reducing agent) to well A-2 while stirring with a clean stirring rod until there is evidence of a reaction.
- 6. Add three drops NaHSO₃ solution to well A-3 while stirring with a stirring rod. Note evidence for a reaction.
- 7. Slowly add NaHSO₃ solution to well A-4, while stirring with a stirring rod, until there is evidence of a reaction.
- 8. MnO₂ is an insoluble brown solid. Mn²⁺ is faintly pink in solution. MnO₄²⁻ is dark green in solution. Based upon these colours, identify the manganese reaction products formed in wells A-2, A-3, and A-4.
- 9. Thoroughly wash your hands before leaving the laboratory.

Data Analysis and Concept Development

- 1. For each well, determine whether the Mn in MnO₄ undergoes oxidation or reduction.
- 2. Fill in the blanks: "MnO₄" acts as _____ (oxidizing or reducing) agent because
- 3. Fill in the blanks: "the other agent in this reaction is _____. It acts as a (oxidizing or reducing) agent because .
- 4. Predict the effect of replacing NaHSO₃ with NaHSO₄. Test your prediction.
- 5. Complete the table involving reactions between KMnO₄ and NaHSO₃ in your lab-book.

	Well#	Change in oxidation state of Mn	Change in oxidation state of S	Electrons gained by Mn	Electrons lost by S
Reaction 1	A-2 (basic)				
Reaction 2	A-3 (neutral)				
Reaction 3	A-4 (acidic)				

- In Reaction 1 of the table above, each MnO₄⁻ takes up one electron, but each HSO₃⁻ gives up two electrons. Account for this difference.
- Predict the stoichiometric ratio of MnO₄ to HSO₃ in Reaction 2 and 3.
- Balance each of the following skeleton reactions based upon the stoichiometric ratios developed in Question 7:

$$OH^{-} + MnO_{4}^{-} + HSO_{3}^{-} ----> MnO_{4}^{2-} + HSO_{4}^{-} + H_{2}O$$

$$H^{+} + MnO_{4}^{-} + HSO_{3}^{-}$$
 ----> $MnO_{2} + HSO_{4}^{-} + H_{2}O$

$$H^{+} + MnO_{4}^{-} + HSO_{3}^{-}$$
 ----> $Mn^{2+} + HSO_{4}^{-} + H_{2}O$

Conclusion

Application

Balance the following redox reactions involved in the production of copper:

1.
$$CuFeS_2 \longrightarrow Cu_2S + FeS + S$$
 (matte)

2.
$$Fe_2O_3 + FeS \longrightarrow FeO + SO_2$$

3. FeO +
$$SiO_2$$
 \longrightarrow FeSiO₃ (Slag formation)

4.
$$Cu_2 S + O_2 \longrightarrow Cu + SO_2$$
 (Blister copper)